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
The screenshot shows the registration form with the following fields and values:

Product family	Mascot Server
Licence type	Permanent
Licensed features	Version 2.6 1 CPUs
Email address*	takaesu@matrixscience.com
Verify email address*	takaesu@matrixscience.com
Full name*	Hiroto Takaesu
Organisation*	Matrix Science Japan
Address line 1*	1-6-10, Higashi Ueno, Taito-ku
Address line 2	
City*	Tokyo
Zip/postcode	110-0015
State/region	
Country*	Japan
Phone number	
Fax number	

A yellow box on the right side of the form contains the text: **\*欄は入力必須** (Required input) and **日本語×、すべて英語で** (No Japanese, all in English).

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of the  
end user, please fill in the end user's details in the  
fields above and your e-mail address here.


Email address (CC)

Matrix Science may send me technical support information by email, such as the availability of new releases, upgrades, and support contracts.

I do not wish to receive any communications from Matrix Science.

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I do not accept the terms of the licence agreement.

I accept the terms of the licence agreement.\*

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Email address (CC)

Matrix Science may send me technical support information by email, such as the availability of new releases, upgrades, and support contracts.
   
 I do not wish to receive any communications from Matrix Science.

**Note:** In either case, Matrix Science will not use your contact information to send unrelated advertising, and will not disclose your details to any third party.

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**MASCOT PROTEIN IDENTIFICATION SYSTEM**

**End-user Licence Agreement**

**IMPORTANT – PLEASE READ CAREFULLY:** This End User Licence Agreement is a legally binding contract between you (either an individual or a single corporate entity) and Matrix Science Limited for the product identified above, which includes computer software, electronic documentation, any printed documentation, and any subsequent updates and supplements (the "Software").

By installing or using the Software, you agree to be bound by the terms of this agreement. If you do not agree to the terms of this agreement, we are unwilling to license the Software to you. In this case, do not install or use the Software. Return the package that included the Software to Matrix Science Limited or their authorised distributor within 30 days of receipt for a full refund.

**1 Licence**


Matrix Science Limited owns the copyright in the Software contained within this package and all other copies which you are authorised by this agreement to make.

This licence is personal to you (either an individual or a single corporate entity) as the purchaser of a licence to use the Software and the licence


I do not accept the terms of the licence agreement.
   
 I accept the terms of the licence agreement.\*

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
The image shows two browser windows. The top window displays the Mascot search engine home page. A yellow callout box highlights the navigation path: Home → Configuration Editor → Database Manager. The bottom window shows the Mascot Configuration page, where 'Database Manager' is circled in red. An orange arrow points from this link to the 'Configuration Editor' link on the home page.

**Mascot Configuration**

<a href="#">Amino Acids</a>	Amino Acid Data
<a href="#">Modifications</a>	Modification definitions
<a href="#">Symbols</a>	Symbols used in chemical formulae
<a href="#">Enzymes</a>	Enzyme definitions
<a href="#">Instruments</a>	Fragmentation Rules
<a href="#">Quantitation</a>	Quantitation Methods
<a href="#">Configuration Options</a>	Global Options in mascot.dat
<a href="#">Database Manager</a>	Sequence databases, Parse Rules and automated downloads

**Mascot Utilities**

- [Database Status](#)  
View the status of all Mascot databases. Links to database statistics, search status, log files.
- [Search Log](#)  
A tabular view of the search log. Can be filtered to find specific search results.
- [Configuration Editor](#)  
Browser-based configuration editors, including security administration and sequence database management.

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The screenshot shows the Mascot Database Manager interface. The main content area is titled 'Databases and spectral libraries'. It contains a table with columns: Name, Mode, Type, Status, and Latest task. The table lists two databases: PRIDE\_Contaminants (predefined, SL, waiting) and SwissProt (predefined, AA, In use). The sidebar on the left has two main sections: 'Database Manager' and 'Library'. Under 'Database Manager', there is a sub-section 'Fasta' with a link 'Enable predefined definition' highlighted in red. Under 'Library', there is a sub-section 'Library' with a link 'Enable predefined definition' highlighted in green. Two yellow callout boxes are overlaid on the screenshot. The first callout box, with red text, points to the 'Fasta' section and says: '配列データベースを追加したい場合、「Fasta」→「Enable predefined definition」をクリック'. The second callout box, with green text, points to the 'Library' section and says: 'ライブラリを追加したい場合、「Library」→「Enable predefined definition」をクリック'.

Name	Mode	Type	Status	Latest task
PRIDE_Contaminants	predefined	SL	(waiting)	Bring online
SwissProt	predefined	AA	In use	Bring online

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Library files will be located in the subdirectory NIST\_Hum directory. The new directory will be created if it does not exist.

Previous Next

**Enable predefined definition**

Library name: NIST\_Human\_IonTrap  
 Sequence directory: C:/inetpub/mascot/sequence  
 Reference database: SwissProt

Please choose a reference database in the spectral library will be taken from the specified Fasta file (the reference database). This will make protein inference more reliable and allows a Protein View report to be displayed for a library hit.

Taxonomy: Homo sapiens (human)

If the selected reference database has taxonomy configured, you can optionally choose a taxonomy for reference accessions. The recommended taxonomy is . . . . . Homo sapiens (human).

Previous Create

**Database: NIST\_Human\_IonTrap**

Copy Delete

Name: NIST\_Human\_IonTrap  
 Using predefined definition: NIST\_Human\_IonTrap  
 Configuration last updated Wed Feb 1 16:03:11 2017

Database type: Spectral library  
 Database directory: C:/inetpub/mascot/sequence/NIST\_Human\_IonTrap/current  
 Filename pattern: NIST\_Human\_IonTrap\_\*.msp

**About to download or copy files**

Database files need to be downloaded or copied before database configuration is completed.  
 Download task is in queue position 1; waiting for downloading to start. See the task queue for more details. (Download task added Thu Feb 2 11:32:14 2017.)

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<対象生物種リスト:DB名XXXXに該当> Environmental , Fungi, Invertebrates, Mammals, Mus, Plants, Prokaryotes, Rodents, Unclassified, Vertebrates

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The screenshot shows the Mascot Database Manager interface. The left sidebar has a red box around 'Running tasks (1)'. An orange arrow points from this box to the 'Running tasks' section of the main content. The 'Running tasks' section shows a table with one task in progress:

Title	Started	Progress
Download files for NIST_Human_HCD <small>dbman_download.pl NIST_Human_HCD</small>	Mon Feb 6 11:07:01 2017	[6.0%] Downloading <ftp://chemdata.nist.gov/download/peptide_library/libraries/human/HCD/2016_05_03/human_hcd_selected.msp.tar.gz> (try 1 out of 5): 915 kB/s (63.20 MB/1.6 GB, 28m 34s remaining)

Below the table is a 'Task log' section with the following data:


Process end time	Command	Last message
[Thu Feb 2 11:41:35 2017]	dbman_download.pl NIST_Human_IonTrap	(success) 'NIST_Human_IonTrap' successfully updated.
[Wed Feb 1 16:18:59 2017]	dbman_download.pl SwissProt	(success) 'SwissProt' successfully updated.

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• Log out  
• Change password  
• Edit settings  
• Current session information

Whether or not Mascot security is enabled, in order to comply with the licensing conditions, access to this Mascot server must be restricted to authorised users. We strongly recommend that your server is protected by a secure firewall. For more information, start from the [WWW Security FAQ](#).

**Mascot Utilities**

Your system administrator can view the status of all the sequence databases. Links to database statuses, search status, log etc.

[Database Status](#)

Search log view of the search log. Can be filtered to find specific search results.

Configuration editors, including security administration and sequence database

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